

## EIA Chlamydia IgM



**EAN Code:** 8595635302879

**Catalog number:** ChM096

**Package size:** 96 tests

**Storage:** 2-8 °C

**Producer:** TestLine Clinical Diagnostics s.r.o.

### Description:

- Microtitre wells are coated with the inactivated and highly purified LPS antigen of strain Chlamydia sp.
- If specific antibodies are present, they bind to the antigen, are labeled by the Conjugate in the following steps and are detected by color reaction with a single component substrate (TMB-Complete).
- The kit allows 96 tests, including controls in a split microtiter plate with color-coded strips and breakable wells.

### Advantages:

- The total assay time is about 2 hours.
- High sensitivity and specificity of the test.
- Kit includes CUT-OFF, Positive Control and Negative Control, Calibrators.
- Semi-quantitative evaluation in the Index of Positivity (IP).
- Ready-to-use, color-coded components.
- Single-component substrate.

- Interchangeable components with the exception of kit specific components (Controls, Conjugate, Plate).

**Application:**

- Dilute samples (1:101).
- Pipette Controls and diluted Samples.
- Incubate at 37°C for 30 min.
- Aspirate and wash the wells 5×
- Pipette Conjugate.
- Incubate at 37°C for 30 min.
- Aspirate and wash the wells 5×
- Pipette Substrate (TMB-Complete).
- Incubate at 37°C for 30 min.
- Pipette Stop Solution.
- Read color intensity at 450 nm.
- Evaluate the test.