


# DS2 & DSX

Fully automated ELISA analyzers



# Quality – flexibility – reliability

The DS2 is a fully open 2-plate analyzer allowing up to 6 parallel tests in 96 wells format according to workflow defined by the diagnostic kit manufacturer.

The DSX 4-plate analyzer is capable of performing up to 12 parallel assays and it is designed to process multiple samples.

Both analyzers consist of one central sample station, pipetting robotic arm, two (DS2) or four (DSX) separate incubators, one absorbance reading station and one washing station.

The analyzers are controlled by an external computer and a software.



## User comfort

- Device installation and training within few hours
- Intuitive software with simple graphical interface
- Time-saving maintenance
- Online help available
- Data export and import in various file formats; quick report and a possibility of sending it by email

## Safety

- Detection system monitoring samples and reagents pipetting
- Detection of liquids volumes, foam, clots or bubbles in a sample
- Washer - liquid-level monitoring
- Acoustic notifications

## Modular system

- Fast repairs, upgrades and device configuration changes
- Automatic recognition and integration of a new module

# Technical specifications DS2 & DSX

Analyzer	DS2	DSX
<b>Capacity</b>	2-plate fully opened	4-plate fully opened
<b>Samples</b>	processing of 96 samples, possibility of two-stage sample dilution in dilution strips	
<b>Primary sample tubes</b>	height 40 – 100 mm diameter 10 – 16 mm	height 40 – 75 mm or 75 – 100 mm diameter 10 – 16 mm
<b>Robotic and pipetting system</b>	1 pipetting arm 2 types of tips (300 µl for samples/controls and 1300 µl for reagents) volume control system	
<b>Tip trays</b>	2 x 108 for samples, 20 for reagents	4 x 108 for samples, 41 for reagents
<b>Pipetting volume</b>	samples 10 – 250 µl, reagents 20 – 1000 µl	samples 10 – 250 µl, reagents 25 – 1000 µl
<b>Pipetting accuracy</b>	< 3% CV (accuracy ± 2%)	
<b>Pipetting arm speed</b>	samples < 15 minutes/plate (50 µl x 96 wells) two-stage dilution < 20 minutes reagents < 2 minutes/plate	samples < 19 minutes/plate (50 µl x 96 wells) two-stage dilution < 26 minutes reagents < 2 minutes/plate
<b>Washing</b>	8-channel wash head, 1 washing station, level sensing	
<b>No. of solutions</b>	2 washing solutions	4 washing solutions
<b>Waste</b>	integrated container (1.5 l)	integrated container (8 l), level sensor
<b>Measurement</b>	1 photometer	
<b>Photometric range</b>	0 – 3,0 OD	
<b>Spectral range</b>	405 – 690 nm	
<b>Incubation</b>	2 heated, independent incubation units	4 heated, independent incubation units
<b>Shaking</b>	continuously or periodically	
<b>Temperature range</b>	RT + 4 °C – 40 °C (±1 °C)	RT + 7 °C – 50 °C (±1 °C)
<b>Hardware</b>	external PC	
<b>Software</b>	operating system Windows XP, bidirectional ASTM interface, XML, etc. archiving patients results, quality control by Levey-Jennings, Westgard rules	
<b>Dimensions (W x H x D)</b>	540 x 660 x 680 mm	1060 x 800 x 910 mm
<b>Certification</b>	CE marked, IVD certified	

## Summary

- Fully open systems (possibility of using all types of diagnostic ELISA kits)
- Possibility of processing Microblot-Array kits
- Bidirectional communication interface
- Easy operation and optimized function settings
- Fast samples and reagents pipetting
- Possibility of two-stage sample dilution, pipetting amount adjustable from 10 to 250 µl
- Optional barcode reader
- Up to 24 or 48 parallel assays per one run

MANUFACTURER:



**DYNEX Technologies Inc.**  
14340 Sullyfield Circle  
Chantilly, VA 20151-1621 USA  
customerservice@dynex.com

DISTRIBUTOR:



**TestLine Clinical Diagnostics Ltd.**

Krizikova 68, 612 00 Brno  
Czech Republic  
[www.testlinecd.com](http://www.testlinecd.com)

**Sales Department:**

+420 549 121 203  
[sales@testlinecd.com](mailto:sales@testlinecd.com)

**Orders:**

+420 541 243 390  
[order@testlinecd.com](mailto:order@testlinecd.com)