Robotikspitzen (im DualPack; Unsteril)

(Made in Germany)

DE - Gebrauchsinformation

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
(Made in Germany)

DE - Gebrauchsinformation

REP UK: MDSS-UK RP Ltd.
6 Wilmslow Road, Rusholme
Manchester M14 5TP
UNITED KINGDOM
E-Mail: info@mdss-ukrp.com

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Robotic tips (in DualPack; non-sterile)

(Made in Germany)

EN – Instructions for use

Intended use of robotic tips from Ritter GmbH:

The robotic tips from Ritter GmbH are used exclusively for the low-contamination dispensing and dosing of liquids in the context of analytical, diagnostic, biotechnological, or microbiological tests. They are intended for single use only and may be used only by qualified and trained specialised/laboratory personnel.

The robotic tips are considered accessories for automated liquid handling platforms, which are themselves marketed as *in vitro* diagnostics (IVD) or general laboratory devices. Together with the necessary reagents or assay kits, these enable the intended use. The instructions for use provided by the manufacturer must be observed with regard to the use, storage, transport, and disposal of the robotic tips.

Non-intended use:

Purification for the purpose of multiple use does not ensure the original functionality of the product and is to be regarded as an unacceptable risk. Re-sterilisation in the case of damaged primary packaging is to be considered an unacceptable risk; the products are to be destroyed.

Delimitation:

STERILE or NON-STERILE, whereby the use of the STERILE form is recommended for diagnostic, biotechnological, or microbiological tests. WITH or WITHOUT FILTER in the inner lumen, whereby the use of the filter variant is recommended for diagnostic, biotechnological, or microbiological tests.

Storage:

- The robotic tips should be stored dry at ambient temperature and protected from direct sunlight.
- Permanent storage in temperatures above 30°C can lead to brittleness of the materials and to loss of performance of the product.
- Storage below ambient temperature can lead to condensation during use. In this case, we recommend storing the products in the laboratory environment for 6–12 h before use.

Application:

When using the products, the following points require your undivided attention:

- Inspection of the information on the carton label. Is the item number correct and has the expiry date not passed?
- Inspection of the shipping carton for damage, moisture, and other signs of improper shipping or storage. If you notice such signs, please refrain from using the tips, and inform Ritter GmbH.
- The primary packaging must be checked for integrity of the tuck flap before opening. If there is any sign of damage to the flap, please refrain from using the tips
- Carefully check the packaging for proper alignment of the rack and the tips before opening. If there are signs of nonconformities, please refrain from using the tips.
- Each rack with tips should be checked for damaged, short, or incorrect tips before use. If there are signs of such defects, please refrain from using the tips.
- For filter tips, check each rack with tips for proper filters: whether there is a filter in each tip and it is not rotated, damaged, or incorrectly positioned. If there are signs of such defects, please refrain from using the tips.
- Each rack with tips must be checked for visible damage before use. If there are signs of damage, please refrain from using the tips.

Disposal:

If there are no other disposal regulations regarding the pipetted liquids that have contaminated the product, the product can be disposed of in the same way as the material described (see product data sheet) according to the nationally applicable legislation.

Manufacturer: Ritter GmbH
Kaufbeurer Strasse 55
D-86830 Schwabmünchen
Tel: +49 8323 5003-45
Email: medical.ritter@avantorsciences.com

REP IN CH: VWR International GmbH
Lerzenstrasse 16/18
CH-8953 Dietikon
Email: gudrun.ruttkowski@avantorsciences.com

REP IN UK: MDSS-UK RP Ltd.
6 Wilmslow Road, Rusholme
Manchester M14 5TP
UNITED KINGDOM
Email: info@mdss-ukrp.com

Robotic tips (in DualPack; non-sterile)

(Made in Germany)

EN – Instructions for use

Meaning of the symbols in use:			
	Catalogue number		Quantity
	The product complies with applicable European requirements		The product is considered an <i>in vitro</i> diagnostic medical device or an accessory thereof
	Sterile product; sterilisation was carried out with ionising radiation		Non-sterile
	Do not use if the packaging is damaged.		Do not reuse
	Protect from sunlight		Store in a dry place
	Global Trade Item Number		Batch number
	Date of manufacture		Use by – date after which the medical device may no longer be used
	Unique product identification		Quantity
	Manufacturer		Internal inspection
	Manufactured in Germany		Observe the instructions for use
	Top (this way up)		Fragile – handle with care

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Tillverkare: Ritter GmbH
REP CH: VWR International GmbH