

CLIA Chlamydia trachomatis IgG

EAN Code: 8595635310201

Catalog number: CL-ChtG100

Package size: 100 tests

Storage: 2-8 °C

Producer: TestLine Clinical Diagnostics s.r.o.



The CLIA immunoassay is based on the principle of magnetic particle separation, and the detection system utilizes flash chemiluminescence with the help of acridinium ester, and the intensity of the light flash is measured in relative light units (RLU).

General information about the kit:

- Each individual test contains all the necessary reagents (ready-to-use) for conducting the examination.
- Calibration is included in the package.
- Results are evaluated using AU/ml.

For more information, please refer to [BiovendorGroup CLIA webpage](#).