

Control set CLIA Borrelia CSF IgG

EAN Code: 8595635312151

Catalog number: CL-BCSFGCON

Package size: 2X20 tests

Storage: 2-8 °C

Producer: TestLine Clinical Diagnostics s.r.o.

- The kit is intended for use together with CLIA kits on the KleeYa[®] analyzer.
- The kit includes positive and negative control sera, which are used to verify the functionality of CLIA tests and the accuracy of obtained results.

For more information, please refer to [BiovendorGroup CLIA webpage](#).