

Control set CLIA Chlamydia pneumoniae IgA

EAN Code: 8595635310836

Catalog number: CL-CPACON

Package size: 2x20 tests

Storage: 2-8 °C

Producer: TestLine Clinical Diagnostics s.r.o.



- The kit is intended for use together with CLIA kits on the KleeYa[®] analyzer.
- The kit includes positive and negative control sera, which are used to verify the functionality of CLIA tests and the accuracy of obtained results.

For more information, please refer to [BiovendorGroup CLIA webpage](#).