


# Microblot- Array

Multiplex diagnostics  
in microtiter plate format

## OVERVIEW




Multiplex testing



ELISA-compatible



Efficiency



Quantitative results



Easy-to-use software




Find out more

# Multiplex diagnostics

Microblot-Array is an immunoblot array in the well of a microtiter plate format designed for efficient multiplex diagnostics. It brings ELISA's high throughput to BLOT testing. The comprehensive evaluation is ensured

by using the TestLine Analytics Software in combination with the Microblot-Array Reader, enabling complex image analysis, including results evaluation and connectivity to LIS.

## So much to discover – in a single well.



## Innovation is part of our portfolio.

<b>Panel</b>	<b>Parameters included</b>	<b>Expected availability</b>
Atypical respirations	Mycoplasma pneumoniae, Bordetella, Chlamydia pneumoniae, Legionella pneumophila	In development
Gastrointestinal diseases	Helicobacter pylori, Yersinia enterocolitica, Autoimmune gastrointestinal diseases	In development
Liver-Kidney diseases	Autoimmune hepatitis, Primary biliary cirrhosis	In development
Autoimmune neurological diseases	Paraneoplastic syndrome, Limbic encephalitis	In development
Herpetic infections	EBV, CMV, HSV, VZV, HHV-6	In development

# 4 steps, countless options.

1

## Unbox


Ready to use kit,  
no dilution needed



2

## Compose

Up to 96 various tests in one plate  
thanks to breakable wells



3

## Process

Automated or manual processing  
with standard ELISA routine



4

## Evaluate

Comprehensive evaluation with  
dedicated MBA Reader and Software





FOLLOW US

BIOVENDOR.GROUP



Diagnostic kits are intended  
for professional use in the laboratory.

# Portfolio. Unique composition, comprehensive diagnostics.

## AUTOIMMUNITY

Code	Products	No. of tests
ANApMA96	Microblot-Array ANA plus*	96

\*Check availability in your country.

## INFECTIOUS SEROLOGY

Code	Products	No. of tests
BpAMA48	Microblot-Array Bordetella IgA	48
BpGMA48	Microblot-Array Bordetella IgG	48
BpMMA48	Microblot-Array Bordetella IgM	48
BGMA096	Microblot-Array Borrelia IgG	96
BMMA096	Microblot-Array Borrelia IgM	96
CoVAMA96	Microblot-Array COVID-19 IgA	96
CoVGMA96	Microblot-Array COVID-19 IgG	96
CoVMMA96	Microblot-Array COVID-19 IgM	96
CMGMA48	Microblot-Array CMV IgG	48
CMMMA48	Microblot-Array CMV IgM	48
EBAMA96	Microblot-Array EBV IgA	96
EBGMA96	Microblot-Array EBV IgG	96
EBMMA96	Microblot-Array EBV IgM	96
HpAMA48	Microblot-Array Helicobacter IgA	48
HpGMA48	Microblot-Array Helicobacter IgG	48
HSGMA48	Microblot-Array HSV 1+2 IgG	48
HSMMA48	Microblot-Array HSV 1+2 IgM	48
CAMA096	Microblot-Array Chlamydia IgA	96
CGMA096	Microblot-Array Chlamydia IgG	96
MyAMA48	Microblot-Array Mycoplasma IgA	48
MyGMA48	Microblot-Array Mycoplasma IgG	48
MyMMA48	Microblot-Array Mycoplasma IgM	48
YAMA048	Microblot-Array Yersinia IgA	48
YGMA048	Microblot-Array Yersinia IgG	48



### TestLine Clinical Diagnostics s.r.o.

Křižíkova 68, 612 00 Brno, Czech Republic

+420 549 121 259 | +420 541 243 390

sales@testlinecd.com

www.testlinecd.com



Company is certified to the quality management system standards ISO 9001 and ISO 13485 for in vitro diagnostics.