

## Gemini

**Catalog number:** 9162800000

**Storage:** Bez teplotního režimu

**Producer:** STRATEC Switzerland AG



### Samples

- Capacity: up to 192 serum tubes
- Flexible: e. g. 144 samples + 8 reagents + 16 controls

### Robotic and pipetting system

- 1 pipetting arm for replaceable tips, 300  $\mu$ l and 1100  $\mu$ l
- Processing capacity: 3 plates
- Pipetting volume: 10-300  $\mu$ l, 301-1000  $\mu$ l
- Volume control system

### Washing

- 1 station
- 3 wash solutions with liquid-level monitoring
- Waste: 10 l container
- Number of channels: 8-channel wash head
- Residual volume: < 2,5  $\mu$ l for U-wells, < 4  $\mu$ l for flat bottom wells

## Measurement

- 1 photometer, 0-3,0 OD, 400-700 nm
- 8 interference filters
- Linearity 0 - 2 OD

## Incubation

- 2 independent incubation units
- Temperature range: RT – 45 °C
- Stability:  $\pm 1,5$  °C
- Possibility of shaking during incubation (different speed setting)

## Software

- Integrated PC
- User friendly
- Bi-lateral LIS connectivity

## Safety

- Liquid-level monitoring
- Advanced monitoring system
- Barcode reader – errors elimination
- Possibility of in-tip dilution

## Dimensions

- Depth: 65 cm; height: 75 cm; width: 125 cm
- Weight: 115 kg

## Certification

- CE marked, IVD certified