

TMB-Complete 2 (1 l)



EAN Code: 8595635303630

Catalog number: TMB0L1

Package size: 1000 ml

Storage: 2-8 °C

Producer: TestLine Clinical Diagnostics s.r.o.

Description:

- Solution contains TMB (3,3', 5,5'-TETRAMETHYLBENZIDINE), hydrogen peroxide (H₂O₂) and stabilizers.
- Solution does not contain any toxic compounds.
- The guaranteed absorbance for 100 µl blank is less than 0.100 at 450 nm.
- Shelf life: 30 months.