

## TMB-Complete 2 (100 ml)



**EAN Code:** 8595635303623

**Catalog number:** TMB100

**Package size:** 100 ml

**Storage:** 2-8 °C

**Producer:** TestLine Clinical Diagnostics s.r.o.

### Description:

- Solution contains TMB (3,3', 5,5'-TETRAMETHYLBENZIDINE), hydrogen peroxide (H<sub>2</sub>O<sub>2</sub>) and stabilizers.
- Solution does not contain any toxic compounds.
- The guaranteed absorbance for 100 µl blank is less than 0.100 at 450 nm.
- Shelf life: 30 months.